# TRUEtrack® Blood Glucose Test Strips

## Instructions for Use

### Intended Use

TRUEtrack Blood Glucose Test Strips are used with TRUEtrack Meters to quantitatively measure whole blood glucose. TRUEtrack can be used for self-testing at home or for professional use. TRUEtrack is for *in vitro* diagnostic use only.

### Test Principle

The TRUEtrack Strip is a plastic Strip containing chemistries and electrodes. TRUEtrack measures glucose by using amperometric technology employing a glucose oxidase reaction. When whole blood or Control is drawn into the tip of a Test Strip, glucose in the sample reacts with chemicals and produces an electrical current. The Meter measures electrical current and calculates amount of glucose. The glucose result is displayed as a plasma value.

### Chemical Composition

Glucose oxidase (Aspergillus sp.) 2.5 units, mediators, buffers and stabilizers.


① **Contact End** End of Strip inserted into Meter.

② **Sample Tip** Pointed end of Strip where sample is drawn into Strip.

Top of Test Strip

②

①

### Caring for Test Strips

- Use Strip quickly after removing from vial. Recap vial right away. Strips left outside of vial too long may give an error message.
- Strips **must** be kept in original vial, with vial cap tightly sealed. NEVER transfer Strips from one vial to another.
- Write date opened on vial label when removing first Strip. Discard all unused Strips in vial after date printed on the Strip vial label, or 120 days after date opened, if either date has passed.
- Store Strip vials in a dry place at room temperature below 86°F/30°C. **DO NOT REFRIGERATE OR FREEZE.** DO NOT store in bathroom or kitchen. DO NOT expose to extreme heat or cold, direct sunlight, or high humidity for any length of time.
- DO NOT reuse Test Strips.
- DO NOT bend, cut or alter Test Strips in any way.

### Important Information

- Test Strips are for *in vitro* testing only. **DO NOT consume.**
- ONLY** use TRUEtrack Test Strips with TRUEtrack Meters and TRUEcontrol Control. Using other Meters or Control may give inaccurate results.
- Do not change your medication or treatment plan based on results from the TRUEtrack Strips and TRUEtrack Meter without advice of Doctor or Diabetes Healthcare Professional.
- Discard used Test Strips and Lancets into an appropriate container (biohazard/for sharp objects).
- Using Test Strips past written or printed dates on the Strip vial label may cause inaccurate results.
- Discard any Strips or vials that appear damaged.
- NEVER** use serum, plasma, or clotted blood when testing with TRUEtrack. **ONLY** use fresh capillary whole blood from the finger or forearm.
- TRUEtrack is not recommended for use with venous samples.**
- When using the forearm for blood sample:<sup>1</sup>
  - Check with your Doctor or Diabetes Healthcare Professional to see if forearm testing is right for you.
  - Results from forearm are not always the same as results from finger.
- Use finger instead of forearm for more accurate results:
  - Within 2 hours of eating, exercise, or taking insulin,
  - If your blood sugar may be rising or falling rapidly or your routine results are often fluctuating,
  - If you are ill or under stress,
  - If your forearm test results do not match how you feel,
  - If your blood sugar may be low or high,
  - If you do not notice symptoms when blood sugar is high or low.

### Quality Control (QC) Testing

There are three quality control tests to let you know that the System is working properly.

An automatic self-test is performed each time a Strip is inserted into a Meter. Upon inserting a Strip into the Test Port, if all segments appear and the Apply Sample Symbol appears in the Display, the Meter is working properly.

#### Quality Control Test: Check Strip

TRUEtrack Check Strip is used to make sure that the Meter is working properly. The Check Strip is a hard plastic strip found in Carrying Case of TRUEtrack Meter. When testing with Check Strip,  is shown in Meter Display if Meter is working properly. See Owner's Booklet for more information on Check Strip testing.

#### Quality Control Test: Control

TRUEcontrol Glucose Control is used to check testing technique and System performance. When Control results fall within ranges found on Test Strip vial label of Test Strips being used, System is working properly and testing technique is good. See Owner's Booklet or TRUEcontrol Glucose Control Instructions for Use for more information on Control Testing.

### Important Information

There are three levels of TRUEcontrol Glucose Control available that contain known amounts of glucose. It is important to perform Control Tests with more than one level of Control to assure your System is working properly and your testing technique is good. For more information on obtaining different levels of Control, call 1-800-803-6025 or 1-954-677-4599.

### Blood Glucose Testing

- Wash area to be lanced. Dry thoroughly.
- Allow Meter and Test Strips to sit at room temperature for 30 minutes before using. If opening vial for the first time, write date opened on vial label.
- Check EXP on Strip vial. Do not use if either date has passed. Discard vial and test with new vial.
- Remove one Strip from vial. Recap vial right away.
- Insert Contact End of Strip, contacts facing up, into Test Port of Meter. Meter turns on. Do not remove Strip from Meter until testing is finished.
- Match Code in Meter Display to Code on Strip vial being used. If Codes do not match, see *Troubleshooting*. Do not use Strips if Codes do not match.
- Obtain blood drop.
- After Apply Sample Symbol appears in Display, place Sample Tip of Strip gently against blood drop and hold. Allow Strip to draw in blood (1.0 µL) until meter beeps. Remove Strip from drop. Dashes appear as Meter is testing. If Meter does not beep 5 seconds after touching Strip to blood drop, see *Troubleshooting*.
- Result is displayed. Record result.
- Remove Strip from Meter. Meter turns off. Discard used Strip and lancet in appropriate container.

### Expected Results for people without diabetes:<sup>2</sup>

	Plasma Blood Glucose Result
Before eating	<100 mg/dL
Two hours after meals	<140 mg/dL
Bedtime	<120 mg/dL

A Doctor or Diabetes Healthcare Professional determines personal target glucose ranges. Always check with your Doctor or Diabetes Healthcare Professional before changing your treatment program.

If you are having symptoms that your glucose is too low or too high, contact your Doctor or Diabetes Healthcare Professional right away.

The TRUEtrack is optimized for testing with fresh capillary whole blood and calibrated to display a plasma result. Plasma results are estimated to be 10%-15% higher than whole blood results.<sup>3</sup>

If comparing results using TRUEtrack Strips to laboratory test results, perform a fingerstick blood test within 15 minutes of the laboratory test. Results are considered accurate when they are within 20% of laboratory results. If you have eaten recently, results using TRUEtrack Strips can be up to 70 mg/dL higher than laboratory results. If glucose tests are performed using capillary blood with the TRUEtrack and venous blood with the laboratory system, capillary blood (finger) may be up to 70 mg/dL (3.9 mmol/L) higher than venous blood (vein), even when each test is performed within 10 minutes after sampling and per manufacturers' instructions.<sup>3</sup>

### Troubleshooting (See Owner's Booklet *Display Messages* for more details.)

If your result is unusually high or low or doesn't match the way you feel, perform a Control Test (see *Quality Control Tests*).

If the Control Test is within range:

- Read *Blood Glucose Testing* again.
- Test your blood glucose again using a new TRUEtrack Strip.

If the results still do not match the way you feel, check with your Doctor or Diabetes Healthcare Professional before changing your treatment plan.

If the results are not within range:

- Check the Expiration Dates. Do not use if past either written date or date printed on Test Strip vial or Control bottle. Test with new Strips/Controls.
- Check for error messages. If an error message appears, follow the Actions in the *Display Message* Section of the Owner's Booklet.
- Check your testing technique. Perform another Control Test.
- Check the temperature. Allow System to reach room temperature 59°-86°F (15°-30°C) before testing.

### Limitations

TRUEtrack products must not be used for testing blood glucose in neonates.

The following **WILL** affect accurate test results:<sup>4</sup>

Abnormally high doses of acetaminophen.

The following **MAY** affect accurate test results:<sup>8,9,10</sup>

Hyperglycemia with hyperosmolarity, with or without ketosis

Thickened blood caused by dehydration

**NOTE:** Similar observations have been reported by other blood glucose system.<sup>5,6</sup>

The following will **NOT** affect accurate test results:<sup>4</sup>

Hematocrit levels between 30% and 55%.

Elevated total blood cholesterol and triglycerides.

Salicylate occurring in expected blood concentrations.

Testing at altitudes up to and including 10,150 feet.

Inaccurate results may occur in severely hypotensive individuals or patients in shock. Inaccurate results may occur for individuals experiencing a hyperglycemic-hyperosmolar state, with or without ketosis. Critically ill patient should not be tested with blood glucose meters.

### Performance Characteristics

#### Accuracy

Accuracy describes how well TRUEtrack results agree with results from an accepted laboratory analyzer. The accuracy of TRUEtrack was assessed against Yellow Springs Instrument Model 23A<sup>4</sup>, (which also uses a glucose oxidase method). Testing was performed by lay users at 4 independent clinical sites. The study demonstrates that TRUEtrack compares well with the laboratory method. A regression analysis of the 229 results is given below:

N	- 229	Slope	- 0.98
Y-intercept	- 5.18	r	- 0.96

#### Precision

Precision describes the variation between TRUEtrack results. Precision studies were performed in a laboratory.

Within Run Precision (whole blood)

	<u>55 mg/dL</u>	<u>82 mg/dL</u>	<u>122 mg/dL</u>	<u>241 mg/dL</u>	<u>500 mg/dL</u>
C.V. (%)	5.1	4.2	3.7	2.7	3.0

Total Precision (Control Solution)

	<u>90 mg/dL</u>	<u>161 mg/dL</u>	<u>286 mg/dL</u>
C.V. (%)	6.3	5.6	4.5

#### Consumer Study

TRUEtrack was tested at diabetes clinics of a large urban area hospital and diabetes care centers. TRUEtrack results were obtained from 212 patients and were compared to parallel results obtained on Yellow Springs Model 2300 STAT. This study shows that patients, when self-testing, and laboratory professionals are able to obtain accurate results when compared to results obtained from a laboratory system.


A clinical study was performed at an urban area hospital to compare forearm and finger stick blood glucose results using the TRUEtrack. TRUEtrack blood glucose tests were performed by patients self-testing (n=100) and on patients' blood by healthcare professionals (n=108). Results demonstrate substantial equivalence between forearm and finger blood glucose values when patients' blood glucose is in steady state and blood glucose tests are performed by patients or healthcare professionals using the TRUEtrack.

**See your TRUEtrack Owner's Booklet for detailed instructions. Call Home Diagnostics, Inc. at 1-800-803-6025 (USA) or 1-954-677-4599 for assistance or visit our website at [www.homediagnostics.com](http://www.homediagnostics.com) for more information.**

**For medical assistance, contact a Doctor or Diabetes Healthcare Professional.**

#### References

- U.S. Food and Drug Administration. *Blood Glucose Meters, Getting the Most Out of Your Meter*. [Electronic Version]. Retrieved June 16, 2009: [www.fda.gov/MedicalDevices/Safety/AlertsandNotices/TipsandArticlesonDeviceSafety/ucm109371.htm](http://www.fda.gov/MedicalDevices/Safety/AlertsandNotices/TipsandArticlesonDeviceSafety/ucm109371.htm).
- Joslin Diabetes Center. *Goals for Blood Glucose Control*. [Electronic Version]. Retrieved June 17, 2009: [www.joslin.org/Beginners\\_guide\\_523.asp](http://www.joslin.org/Beginners_guide_523.asp).
- Sacks, DB.: "Carbohydrates" Burtis, C.A. and Ashwood, E.R. (eds.), Tietz Textbook of Clinical Chemistry. Philadelphia: W.B. Saunders Company (1994), 959.
- Data on file.
- Young, D.S. *Effects of Drugs of Clinical Laboratory Tests*, Third Edition. Washington, D.C.: AACC Oress (1990), pp. 3-168, 3-180.
- Sacks, D.B., "Carbohydrates," Burtis, C.A., and Ashwood, E.R. (ed.), Tietz Textbook of Clinical Chemistry. Philadelphia: W.B. Saunders Company (1994), 959.
- Yellow Springs Instrument Model 23A Glucose Analyzer Instruction Manual, "Principles of Operation", pp. 1-2 and "Performance Specifications", pp. 3, 1986.
- Atkins, S.H.; Dasmahapatra, A.; Jaker, M.A.; Chorost, M.I.; Redd, S., Fingerstick Glucose Determination in Shock. *Annals of Internal Medicine*, 114:1020-1024 (1991).
- Sander, M.; Low-Beer, T., Misleading Capillary Glucose Measurements, *Practical Diabetes*, 7:210 (1990).
- Wickham, N.W.R.; Achar, K.N.; Cove, D.H., Unreliability of Capillary Blood Glucose in Peripheral Vascular Disease, *Practical Diabetes*, 3:100 (1986).

A3HD03R12

TRUEtrack

TRUEtrack

Manufactured by:

HOMEdiagnostics™

Fort Lauderdale, FL 33309 U.S.A.  
www.homediagnostics.com

TRUEtrack®

© 2009 Home Diagnostics, Inc.  
TRUEtrack and the Home Diagnostics logo are trademarks of Home Diagnostics Inc.  
A3HD03 Rev. 12

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### Referencias

- U.S. Food and Drug Administration. *Blood Glucose Meters, Getting the Most Out of Your Meter*. [Electronic Version]. Retrieved June 16, 2009: [www.fda.gov/MedicalDevices/Safety/AlertsandNotices/TipsandArticlesonDeviceSafety/ucm109371.htm](http://www.fda.gov/MedicalDevices/Safety/AlertsandNotices/TipsandArticlesonDeviceSafety/ucm109371.htm).
- Joslin Diabetes Center. *Goals for Blood Glucose Control*. [Electronic Version]. Retrieved June 17, 2009: [www.joslin.org/Beginners\\_guide\\_523.asp](http://www.joslin.org/Beginners_guide_523.asp).
- Sacks, DB.: "Carbohydrates" Burtis, C.A. and Ashwood, E.R. (eds.), Tietz Textbook of Clinical Chemistry. Philadelphia: W.B. Saunders Company (1994), 959.
- Data on file.
- Young, D.S. *Effects of Drugs of Clinical Laboratory Tests*, Third Edition. Washington, D.C.: AACC Oress (1990), pp. 3-168, 3-180.
- Sacks, D.B., "Carbohydrates," Burtis, C.A., and Ashwood, E.R. (ed.), Tietz Textbook of Clinical Chemistry. Philadelphia: W.B. Saunders Company (1994), 959.
- Yellow Springs Instrument Model 23A Glucose Analyzer Instruction Manual, "Principles of Operation", pp. 1-2 and "Performance Specifications", pp. 3, 1986.
- Atkins, S.H.;; Dasmahapatra, A.; Jaker, M.A.; Chorost, M.I.; Redd, S., Fingerstick Glucose Determination in Shock. *Annals of Internal Medicine*, 114:1020-1024 (1991).
- Sander, M.: Low-Beer, T., Misleading Capillary Glucose Measurements, *Practical Diabetes*, 7:210 (1990).
- Wickham, N.W.R.; Achar, K.N.; Cove, D.H., Unreliability of Capillary Blood Glucose in Peripheral Vascular Disease, *Practical Diabetes*, 3:100 (1986).

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HOMEdiagnostics™

Fort Lauderdale, FL 33309 U.S.A.  
[www.homediagnostics.com](http://www.homediagnostics.com)

TRUEtrack®